UROLOGIX[®]

PROSTIVA[®] RF THERAPY

0° Telescope 8099

15° Telescope 8099TU15

Télescope à 0° 8099

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0°-Teleskop 8099

15°-Teleskop 8099TU15

Endoscopio a 0° Modelo 8099

Endoscopio a 15°

Modelo 8099TU15

0°-telescoop 8099

15°-telescoop 8099TU15

**Ottica endoscopica a 0°
modello 8099**

**Ottica endoscopica a 15°
modello 8099TU15**

0° teleskooppi 8099

15° teleskooppi 8099TU15

0°-behandlingsoptik 8099

**15°-behandlingsoptik
8099TU15**

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Prostiva RF Therapy Telescope

Introduction

The instructions below describe the general procedure for safe handling of the Prostiva RF Therapy Telescopes. The telescope is a reusable device. Before you use the telescope and for each subsequent use, it must be cleaned, decontaminated, and sterilized. These instructions are designed to assist you in cleaning, sterilizing, caring for, and storing your telescopes. Please read them carefully. By following these instructions, your telescopes will provide you with a long service life.

Indications for Use

The telescope is used to position the Model 8929 Prostiva RF Hand Piece through the urethra into the prostate. This telescope shall only be used in conjunction with the Prostiva Hand Piece.

Safety Information

For a list of contraindications, warnings, precautions and adverse events refer to the Prostiva RF Therapy System User Guide.

Instructions for Use

Inspection

Upon receipt, immediately inspect the telescope for signs of damage. Also inspect the Prostiva RF Therapy Telescope before and after every procedure. Dents along the shaft and scratches on the lens are the most noticeable signs of damage. Look through and rotate the telescope to check the clarity of the view. Scratches, fingerprints and debris will cause an impaired view. Additionally, debris may act as an abrasive and scratch the lens. To make certain a telescope is in good working order, observe the reflected light on the surfaces of the ocular and objective lenses. These surfaces should be smooth and shiny.

Inspect the integrity of the fiber optics by holding the light post towards a light and observing the distal tip. The fiber optics are intact when the light is evenly distributed. Darkened areas indicate broken fibers. Foggy images may result from moisture entering a damaged seal around the lens. A partially or completely obstructed view may be the result of damaged lenses within the

telescope shaft or eyepiece. Do not use a telescope that shows signs of damage.

If it is necessary to dispose of the telescope, discard in accordance with local environmental regulations.

For assembly, installation and connection instructions, refer to the Prostiva RF Therapy System User Guide.

Handling

The Prostiva RF Therapy Telescope consists of an eyepiece with glass ocular lens, a fiber optic light cable connection with screw-on type cable connectors to accept several manufacturers' light cables, and a stainless steel shaft containing the rod lens system and objective lens.

As with all optical instruments, the telescope is a very delicate piece of endoscopic equipment. To ensure the long-lasting quality of the telescope, please follow the handling instruction described below.

- To prevent breakage, the telescope should be supported by firmly grasping the eyepiece end. Never handle by the distal end alone.
- Never bend the stainless steel shaft. This could lead to breaks or cracks in the rod lens system.
- Handle the telescope with care. Hard knocks, particularly at the distal end, may result in damage or cracks in the telescope and allow liquid, steam and other material to penetrate. If so, the image will appear foggy or have an unclear visual area.
- To prevent damage when cleaning, disinfecting and sterilizing the telescope, it is recommended that they be handled separately from other instruments.

Cleaning

General Considerations – Adequate reprocessing is contingent upon the thoroughness of cleaning. Cleaning is the initial step and sterilization occurs later in reprocessing and is intended to kill microorganisms to reduce the likelihood of transmission and possibilities of infection. To ensure acceptable reprocessing, there should be no delay between the steps in this document.

Bloodborne Pathogens – Universal precautions for handling this device after use should be observed by all hospital personnel according to OSHA Standard 29 CFR 1910.1030, Occupational Exposure to Bloodborne Pathogens. Protective gloves and eyewear should be worn throughout the entire cleaning process.

Thorough Cleaning of this Device – It is critical that this device is thoroughly cleaned before its first use, and before all subsequent uses, following the specific cleaning steps listed in this document. Ineffective or incomplete cleaning can hinder subsequent sterilization activities.

Automated Cleaning –

 Caution: Prostiva RF Therapy Telescopes **MUST NOT** be cleaned in an ultrasonic bath or automatic washer. These cleaning methods can damage the telescope. This can make it nonfunctional and require that it be discarded.

Cleaning Agents – The use of a neutral pH enzymatic cleaning agent is recommended. This type of cleaner has been shown to effectively remove organic soils, such as blood and urine, from this device. The following cleaning instructions have been validated using Enzol enzymatic detergent solution. In addition, the following equipment (or equivalent) is recommended:

- Lumen brush (2), soft nylon bristled, 3.0-6.0 mm (1/8-1/4") diameter (for cleaning adaptors)
- Instrument brush (2), soft nylon bristled, 18.0-20.0 cm (7-8") length (for cleaning scopes)
- Syringe, needleless (60.0 mL or larger)

The following cleaning agents, solutions, or tools should **NOT** be used when cleaning this device:

- Saline solution
- Alkaline cleaning agents
- Solutions containing chlorine (eg, bleach) or aldehydes (eg, gluteraldehyde)
- Formalin, mercury, chlorides, bromides, iodides, or ringers solution
- Metal brushes or scouring pads

Cleaning and Rinsing Water – Use warm tap water to clean this device.

Deionized water should be used for the final rinse step to prevent mineral deposits on the device surfaces.

Manual Cleaning Instructions

Cleaning Instructions: Point of Use – DO NOT allow this device to dry after use and prior to cleaning.

To prevent blood, urine and other contaminants from drying on the telescope, it should be placed in a container with water immediately after use. Do not put other instruments on top of the telescope.

1. Disassemble any light cable adaptors being used from the telescope.
2. Rinse the telescope and light cable adaptors for a minimum of 3 minutes with a steady stream of warm tap water while thoroughly scrubbing with a soft, nylon bristled brush to remove gross contaminants.
3. Using warm tap water (>18 °C >70 °F), prepare an enzymatic cleaning solution according to the detergent manufacturer's instructions (Enzol=1.0 oz/ gallon).
4. Place telescope and light cable adaptors in the cleaning solution, completely submerged, and soak for a minimum of 2 minutes.
5. After soaking, and while still submerged, thoroughly scrub with a second, soft, nylon bristled brush to gently remove all remaining visible contamination from the telescope and light cable adaptors. Flush challenging areas (including the connections near the ocular lens assembly and adaptors) using a needleless syringe completely filled with warm tap water. Repeat scrubbing and flushing as necessary.
6. After cleaning, thoroughly rinse the telescope and light cable adaptors with de-ionized water (>13 °C >55 °F) for a minimum of 2 minutes to remove detergent residue.
7. Use a clean, absorbent, low-lint wipe to dry the telescope and light cable adaptors.
8. Visually inspect the telescope and light cable adaptors:
 - for signs of damage, such as dents along the shaft or scratches on the lenses. Do

not use a telescope that shows signs of damage.

- for any visible contamination. Do not proceed with sterilization of a contaminated device. If any contamination is still present, repeat the above cleaning steps.

9. Upon completion of cleaning, check the telescope image quality for darkened, obstructed, or foggy views. If there are no signs of damage, proceed with sterilization of the telescope.

Sterilization

Due to variations in hospital sterilizers, precise instructions for sterilization or aeration cannot be given here. If further information is necessary regarding the procedure to be used, contact the manufacturer of the sterilizer unit. Use biological indicators or other acceptable methods to assist in validating the effectiveness of the hospital's sterilizer unit.

 **Caution:** Do not attempt to sterilize the telescope before its first use or after subsequent uses without first thoroughly cleaning and removing any organic material, blood or cleaning solution. Failure to properly prepare the telescope for sterilization will adversely affect the sterilization process.

Routine autoclave or ethylene oxide sterilization is recommended for sterilization of the telescope.

The telescope may be sterilized by the following methods:

Autoclave Sterilization

The telescope can be steam autoclaved in accordance with the following instructions:

 **Caution:** Never autoclave the telescope without wrapping with appropriate covering. Autoclaving the telescope without appropriate covering (eg, flash autoclaving in the OR) can damage the telescope, which can make it non-functional and require that it be discarded.

 **Caution:** Do not immediately expose the telescope to air after removal from the autoclave. Never attempt to cool the telescope by pouring cool liquid over the scope. Sudden changes in temperature may fracture

the glass components of the telescope and require that it be discarded.

 **Caution:** Do not allow the telescope to come into contact with any metal during steam autoclaving. Contact with metal can adversely affect the telescope during steam autoclaving, which can make it non-functional and require that it be discarded.

- Clean and dry the telescope by following the Cleaning Instructions on page 12.
- Carefully place the telescope and light cable adapters into the sterilization container(s).
- Double wrap the container with the appropriate wrapping.
- Autoclave at a maximum temperature of 132 °C (270 °F) for 15 minutes at 223 kPa (32 psi).
- When the autoclave cycle is finished, remove the container(s) from the autoclave and allow the telescope and adaptors to cool to room temperature before removing the top of the container(s).

Ethylene Oxide (EtO) Gas Sterilization

- Clean and dry the telescope following the instructions given in the Cleaning Instructions on page 12.
- Telescopes may be aerated without damage for any desired time period at temperatures up to 68.3 °C (155 °F). Follow sterilizer manufacturer's instructions regarding exact temperatures for your aerator.

Temperature	53 °C (127.4 °F)
Vacuum	533 mm Hg (21 in Hg)
Pressure	96.5 kPa (14 psi)
Humidity	70%
Pre-conditioning	1 hour
Carrier	HCFC
Concentration	600 Mg/L
Exposure	2 hours
Aeration	55 °C /12 hours (131 °F/12 hours)

 **Caution:** During gas sterilization, the telescopes should not come into direct contact with metal. Contact with

metal can adversely affect the telescope during gas sterilization, which can make it non-functional and require that it be discarded.

Sterrad®

The Prostiva RF Therapy Telescope may be sterilized by the Sterrad 100 sterilization system. Place the telescope in the Sterrad 100 sterilization system. Place the telescope in the Sterrad instrument tray and wrap with polypropylene sterilization wrap. Follow directions provided with the Sterrad system.

Telescopes sterilized in the Sterrad 100 system may exhibit cosmetic changes; these changes do not affect the functionality of the device.

Chemical Disinfection

The telescope may be chemically disinfected using a high level disinfecting solution which is specially made for endoscopes (eg, Cidex).

Caution: Only EPA registered hard surface disinfection solutions capable of high level disinfection should be used. Solutions with a strong acid or alkaline pH are not recommended. Follow solution manufacturer's recommendations for proper dilution and recommended minimum exposure times for the Telescopes. **DO NOT EXCEED RECOMMENDED MAXIMUM EXPOSURE TIME OR 45 MINUTES WHICHEVER IS LESS.** Prolonged exposure may result in tarnishing or corrosion.

- A plastic basin should be used to avoid scratching the telescope and to eliminate electrolytic corrosion which can occur when dissimilar metals are soaked in the same solution.
- Clean and dry the telescope following the instructions given in the Cleaning Instructions on page 12.
- For chemical disinfection, carefully follow solution manufacturer's instructions. Do not allow the telescope to remain for more than 45 minutes in any solution, including sterile water.
- After removal from the disinfection solution, rinse thoroughly with sterile water.

- Dry the telescope with a sterile cloth.
- Apply alcohol to the fiber surfaces to complete the drying process.

 **Caution:** Alcohol must be applied carefully, paying special attention to the fiber surfaces in the light piece to complete the drying process. Any residues of disinfectant or cleaning agents in the light inlet post could burn into it when the fiber optic light cable is connected and this would considerably impair light transmission.

Important Notice – Equipment Limited Warranty (U.S. Customers only)

- A. This Limited Warranty provides the following assurance to the purchaser of the Prostiva RF Therapy Telescope (8099 or 8099TU15), hereafter referred to as “Equipment”:
- (1) Should the Equipment fail to function within normal tolerances due to a defect in materials or workmanship within a period of 90 days, Urologix will, at its option:
 - (a) repair or replace any part or parts of the Equipment;
 - (b) issue a credit to the purchaser equal to the Purchase Price, as defined in Subsection A(2), against the purchase of the replacement Equipment or
 - (c) provide a functionally comparable replacement Equipment at no charge.
 - (2) As used herein, Purchase Price shall mean the lesser of the net invoiced price of the original, or current functionally comparable, or replacement Equipment.
- B. To qualify for Limited Warranty set forth in Section A(1), these conditions must be met:
- (1) The Equipment must be returned to Urologix within thirty (30) days after discovery of the defect. (Urologix may, at its option, repair the Equipment on site).
 - (2) The Equipment must not have been repaired or altered outside of Urologix’ factory in any way which, in the judgment of Urologix, affects its stability and reliability. The Equipment must not have been subjected to misuse, abuse or accident.
- C. This Limited Warranty is limited to its express terms.

In particular:

- (1) Except as expressly provided by this Limited Warranty, Urologix is not responsible for any direct, incidental or consequential damages based on any defect, failure or malfunction of the EQUIPMENT, whether the claim is based on warranty, contract, tort or otherwise.
- (2) This Limited Warranty is made only to the purchaser who uses the Equipment. AS TO ALL OTHERS, Urologix makes no warranty, expressed or implied, including but not limited to, any implied warranty of merchantability, or fitness for a particular purpose, whether arising from statute, common law, customer or otherwise. No express or implied warranty to the patient shall extend beyond the period specified in A(1) above. This Limited Warranty shall be the exclusive remedy available to any person.
- (3) The exclusions and limitation set out above are not intended to, and should not be constructed so as to contravene mandatory provisions of applicable law. If any part or term of this Limited Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of the Limited Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Limited Warranty did not contain the particular part or term held to be invalid. This Limited Warranty gives the purchaser specific legal rights. The purchaser may also have other rights which vary from state to state.
- (4) No person has any authority to bind Urologix to any representation, condition or warranty except this Limited Warranty.

This Limited Warranty is provided by Urologix, Inc., 14405 Twenty First Avenue North, Minneapolis, MN 55447-4699. It applies only in the United States. Areas outside the United States should contact their local distributor for exact terms of the Limited Warranty.

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Urologix and Prostiva are registered trademarks of Urologix, Inc.
Sterrad and Enzol are registered trademarks of Johnson & Johnson Corp.
All Rights Reserved
© Urologix, Inc. 2012